

# Dishes

## CS Receiver

- Multi-purpose - can be used for various functions such as carrying instruments and dressings.
- Larger capacity than traditional kidney dish.
- Steady when standing due to flat bottom surface.
- Strong sides so able to support weight of contents when full.
- Easy to store (product is stackable) and separate.
- Easy to carry and hold.
- Ergonomically shaped to accommodate the patient's neck.



### Specification

Spec No: 67442  
Capacity: 900ml  
Material: Paper Pulp  
Product Size: 260L x 160W x 40H mm  
Weight: 30g ± 4g  
Fluid Resistance: British Standard: PAS29:1999  
Fluid: Water  
Resistance Time: >4 Hours  
Resistance Temp: 35°C ± 3°C  
Disposal: By maceration in accordance with PAS29:1999

### Packaging

Units Per Carton: 300  
Carton Weight: 9.57kg (average weight)  
Carton Dimensions: 550L x 336W x 595H mm  
Cartons Per Pallet: 12  
Pallet Size: 1200L x 1000W x 1260H mm (approx.)  
Box Weight: 0.57kg  
Pallet Weight: 123kg  
Carton Volume: 0.11m<sup>2</sup>  
20ft Container: 273 Cartons  
40ft Container: 546 Cartons